Instructions for use
TSK ACECUT BIOPSY NEEDLE
TSK ACECUT BIOPSY NEEDLE with Co-axial Guide Needle

The instructions for use covers TSK ACECUT BIOPSY NEEDLE (hereinafter referred to as "the device") and TSK ACECUT BIOPSY NEEDLE with Co-axial Guide Needle [with TSK Co-axial Guide Needle (hereinafter referred to as "introducer needle") attached].

1 Intended purpose

The device is a biopsy needle used to puncture the human body and collect tissue.

2 Indication

The device is indicated for use when a biopsy sample is required for inspection, treatment, or diagnosis.

3 Intended users

The device is intended for use by doctors or medical professionals skilled in the procedure, and not intended to be used by others. Do not use the device unless you are a doctor or a medical professional who is skilled in the procedure.

4 Target patients

The device is intended for use in patients 12 years of age and older of any gender, excluding cases of contraindications. The device should not be used in pregnant and lactating women.

5 Performance characteristics

- Punctures the human body
- Collects tissue for inspection, treatment or diagnosis

6 Clinical benefits

The biopsy needle is a simple ancillary device which is widely used in many clinical procedures; the clinical benefit of the device is to have a positive impact on patient management:

- Allow management of patients' medical conditions by providing a means of collecting tissue samples locally in order to meet patients' medical needs
- Allow medical professionals to manage medical conditions that require removal of tissue samples, and obtain tissue samples for diagnostic purposes

7 Undesirable side effects

Possible adverse events associated with the use of the device are listed in the following:

- pain, post-procedural bleeding, hematoma, vasovagal reactions, subcutaneous emphysema, arteriovenous fistula due to puncture, edema, bruising, pseudoaneurysm formation, infection, inflammation, sepsis, fever, perforation of unintended tissue, needle stick injuries, breaking of needle or needle tip damage, device self-activation or other malfunction, failure to obtain sample, tumor seeding along needle tract
- side effects associated with local anesthesia (in rare cases, hypersensitivity reactions to local anesthetics such as hypotension, clouding of consciousness (loss of consciousness), and convulsions may occur

Additionally, the following site-specific adverse events could occur:

- Prostate biopsy:
 - Hematuria
 - bloody stool/serous fluid
 - dysuria
- Kidney biopsy:
 - Hemorrhage
 - bacterial infection
 - arteriovenous fistula
- Liver biopsy:
 - fever, abdominal pain (caused by skin infection, bile leakage due to puncture of the biliary tract, cholangitis, etc.)
 - transient liver dysfunction
 - pneumothorax (a hole in the lung during specimen collection may cause a pneumothorax)
 - in the case of biopsy of liver tumor in the liver, spreading of cancer cells (rare case)
- Lung biopsy:
 - air embolism due to puncture operations in the lungs, leading to cerebral infarction and myocardial ischemia
- Thyroid biopsy:
 - Dysphagia
- Breast biopsy:
 - darkening of the surrounding area due to posterior bleeding
 - pain

8 Contraindications

- patients who are receiving anticoagulant therapy, or have bleeding disorders or bleeding tendency

9 Warning

9.1 General precautions

1.	Do not reuse. [Otherwise, there is a risk of infection.]
2.	Do not resterilize. [Re-sterilization does not assure device sterility and could damage or degrade the device.]
3.	Do not use the device under MRI. [Otherwise, there is a risk of MRI accidents.]
4.	Before using the device, read and follow the instructions for use.
5.	If the patient or the device shows any abnormality, discontinue use immediately and take necessary measures.

6.	Check the device's appearance before use to ensure no damage (scratches, cracks, deformation, discoloration, and others).
7.	Refer to the instructions if you use the device with pharmaceuticals or other medical devices.
8.	Patient management after biopsy depends on the biopsy technique and the physiological state of the individual patient. Monitor the patient's vital signs and strive to prevent and treat adverse events associated with the biopsy procedure.
9.	Do not use the device for any purpose other than those described in the "instructions for use."
10.	Do not modify the device.
11.	Operation practicing results in damage to the needlepoint and poor sharpness. Use samples for operation practice.
12.	The device could contain the following substances, in a concentration above 0,1 wt%. <ul style="list-style-type: none"> • Cobalt: CAS No. 7440-48-4, EC No. 231-158-0 Current scientific evidence confirms that medical devices manufactured from cobalt alloy or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

9.2 Precautions for use

- Open and remove with the appropriate aseptic technique. [Otherwise, there is a risk of infection.]
- Use the device immediately after opening the package. [Otherwise, bacteria could multiply and cause infection.]
- Do not touch the needle directly with your hands. [Otherwise, there is a risk of infection.]
- Confirm the selection of gauge size and stroke length in advance. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- When using an introducer needle, use an introducer needle that is compatible with the device. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- When using a biopsy needle guide such as an ultrasonic probe in combination, do not forcibly set the needle on a guide whose shape does not fit it. [Otherwise, there is a risk of device damage.]
- The amount of tissue sampling depends on the nature of the target tissue.
- Be careful when removing or attaching the protector. [Otherwise, there is a risk of needle stick injury.]
- The location of the puncture site should be confirmed by diagnostic equipment and palpation.
- Check the puncture route carefully when inserting the needle. [Otherwise, there is a risk of damaging unintended areas.]
- Do not deflect the needle during operation. [Otherwise, there is a risk of device damage.]
- Be careful not to let the needlepoint touch anywhere other than the intended contact surface. [Otherwise, there is a risk of the needlepoint damage.]
- There have been reports of air embolism due to puncture operations in the lungs, leading to cerebral infarction and myocardial ischemia. Thus, if the patient develops any of these symptoms, keep the patient's head in a low position, diagnose with CT and others, and take appropriate measures immediately. In severe cases, consider hyperbaric oxygen therapy without delay.
- Do not use the device to inject or drain fluids into or out of the body.
- After the procedure, take appropriate measures such as hemostasis.
- When wiping off the inner or outer needle, be careful not to get needle stick injuries and do not overload the device. Also, ensure that no foreign substance adheres to the needle parts. [Otherwise, there is a risk of device damage.]
- When performing repeated biopsies, check the appearance and ensure no damage (scratches, cracks, deformation, discoloration, and others) before use.
- Note that after biopsy, the blood and disinfectant (iodine tincture and others) could clot, and the levers could become stuck.


9.3 Other precautions

- The user must report all adverse events resulting from the use of the device to the manufacturer and to the competent authorities of the member state where the user or patient is.
- After use, dispose of the device safely and with care to prevent the risk of infection, following institutional, local, regional and national guidelines.

10 Device specifications


10.1 Configuration (typical diagram)

TSK ACECUT BIOPSY NEEDLE

TSK Co-axial Guide Needle

Only model of TSK ACECUT BIOPSY NEEDLE with Co-axial Guide Needle is attached.

10.2 Material

Raw material for parts that contact the blood, body fluids, mucous membranes, and others is below.

- Inner and outer needles: stainless steel

The device is ethylene oxide sterilized.

11 Instructions for use

11.1 Procedure of using introducer needle

If not use the introducer needle, perform from 4).


- 1) Remove the protector of introducer needle.
- 2) Puncture the introducer needle towards the marginal zone of the target site. If there is resistance and the needle does not penetrate as expected during insertion, do not forcibly insert the needle. [Otherwise, there is a risk of device damage.]
- 3) Remove the inner needle from the outer needle of the introducer needle. When removing the inner needle, hold the outer needle hub and turn the inner needle hub to unlock it.

11.2 Preparation for biopsy


- 4) Remove the protector from the device.
- 5) Set the safety switch to position "lock" with your thumb and index finger.

- 6) Confirm that the sampling lever is in the closed position. The closed position means that the sampling lever is in the state shown in the figure below.

- 7) Pull up the cover from the body until it "clicks" into place. If the sampling lever also moves, the lever is not in the closed position, so stop pulling up the cover and check that the lever is in the closed position again. Do not pull up while the lever is not in the closed position. [Otherwise, there is a risk of device damage.]

- 8) After pulling up the cover, replace it with its original condition before performing subsequent operations. [Otherwise, there is a risk of device damage.]


11.3 Biopsy procedure

If not use the introducer needle, perform from 10).

- 9) Insert the needle of the device into the outer needle of the introducer needle.
- 10) Puncture the device towards the marginal zone of the biopsy site. If there is resistance and the needle does not penetrate as expected during insertion, do not forcibly insert the needle. [Otherwise, there is a risk of device damage.]

Fully automatic firing (one-step sampling)

- 11) Set the safety switch to position "fire."


- 12) Push the Trigger button all the way in to advance the inner/outer needle. Push sufficiently. [Otherwise, they could fail to move forward correctly.]
- 13) Continued in 15)

Semi-automatic firing (two-step sampling)


- 11) Set the safety switch to the line between fire and lock ("—" position).
- 12) Push the Trigger button to advance the inner needle.
- 13) Set the safety switch to position "fire."
- 14) Push the Trigger button to advance the outer needle.

11.4 Removal of tissue

- 15) After removing the needle, pull the sampling lever down to the open position, as shown in the figure below, and remove the tissue from the tissue sampling groove.

- 16) After removal, return the sampling lever securely to the closed position.

- 17) If you repeat the biopsy, return to 11.2 Preparation for biopsy.

11.5 Removal the introducer needle

If not use the introducer needle, the following procedure is not performed.

- 18) Insert the inner needle into the outer needle of the introducer needle. Do not remove the outer needle alone. [Otherwise, there is a risk of device damage.]
- 19) While holding the outer needle hub, turn and lock the inner needle hub.
- 20) Remove the introducer needle.

12 Storage

- Keep dry. Keep away from direct sunlight, high temperature, and humidity. Store in a clean place.

13 Expiry date

- For the useful life of the device that has been stored properly as per the storage instructions, refer to the use-by date printed on each individual package.
- Store the device properly and do not use after the expiry date.

14 Manufacturer's information

TSK Laboratory, Japan
 1510-1 Soja-Machi
 Tochigi-Shi, Tochigi-Ken
 328-0002 JAPAN
 TEL: +81-282-27-8426
 FAX: +81-282-25-6511
 URL: www.tsklab.co.jp
 E-mail: tsk@tsklab.co.jp

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**TSK ACECUT BIOPSY NEEDLE
TSK ACECUT BIOPSY NEEDLE with Co-axial Guide Needle**

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