Instructions for use

TSK BMB SURELOCK TSK BMB

1 Intended purpose

The TSK BMB SURELOCK is a needle biopsy system for taking histological and cytological bone marrow tissue specimens for inspection, diagnosis, or treatment of patients.

2 Indications

The device is a needle biopsy system for taking histological and cytological bone marrow tissue specimens.

3 Intended users

The device is intended for use by medical doctors and healthcare professionals in healthcare facilities.

4 Target patients

The device is intended for use in human of any gender, excluding cases of contraindications. The device should not be used in pregnant and lactating women.

5 Performance-characteristics

The device is needle biopsy systems to collect histological and cytological bone marrow tissue samples for the examination, diagnosis, or treatment of patients.

6 Clinical benefits

- To provide a means of locally harvesting cells to manage a patient's health condition and meet the patient's medical needs.
- To enable healthcare workers to collect cells for diagnostic purposes and manage required health conditions.

7 Undesirable side effects

Possible adverse events associated with the use of the device are listed in the following:

- organ damage
- injury
- injury due to product breakage
- bacterial infection/fever/inflammation/secondary infection/morbidity
- bleeding/hemorrhage
- cardiac tamponade followed by death

8 Contraindications

- Do not use for patients who are receiving anticoagulant therapy, or have bleeding disorders or bleeding tendency.
- Do not use for the heart, the central circulatory system and the central nervous system.

9 Prohibitions

- Do not reuse. [Otherwise, there is a risk of infection.]
- Do not resterilize. [Re-sterilization does not assure device sterility and could damage or degrade the device.]

10 Warning

10.1 General precautions

1.	Do not use the device under MRI. [Otherwise, there is a risk of MRI accidents.]
2.	Before using the device, read and follow the instructions for use.
3.	Do not use the device unless you are a medical doctor or a health care professional skilled in the procedure.
4.	If the patient or the device shows any abnormality, discontinue use immediately and take necessary measures.
5.	Check the device's appearance before use to ensure no damage (scratches, cracks, deformation, discoloration, and others).
6.	Do not use the device if the sterile bag damage (scratched, cracks, deformation, discoloration, and others) or the unintended opening before use are shown.
7.	Refer to the instructions if you use the device with pharmaceuticals or other medical devices.
8.	Patient management after biopsy depends on the biopsy technique and the physiological state of the individual patient. Monitor the patient's vital signs and strive to prevent and treat adverse events associated with the biopsy procedure.
9.	Do not use the device for any purpose other than those described in the "instructions for use."
10.	Do not modify the device.
11.	Operation practicing results in damage to the needlepoint and poor sharpness. Use samples for operation practice.
12.	Use a syringe that comply with ISO 80369-7.

10.2 Precautions for operation

- Open and remove the device with the appropriate aseptic technique. [Otherwise, there is a risk of infection.]
- Use the device immediately after opening the package. [Otherwise, bacteria could multiply and cause infection.]
- Do not touch the needle directly with your hands. [Otherwise, there is a risk

of infection.]

- Confirm the selection of gauge size and stroke length in advance. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- When using a biopsy needle guide such as an ultrasonic probe in combination, do not forcibly set the needle on a guide whose shape does not fit it. [Otherwise, there is a risk of device damage.]
- The amount of tissue sampling depends on the nature of the target tissue.
- Be careful when removing or attaching the protector. [Otherwise, there is a risk of needle stick injury.]
- The location of the puncture site should be confirmed by diagnostic equipment and palpation.
- Check the puncture route carefully when inserting the needle. [Otherwise, there is a risk of damaging unintended areas.]
- Perform each operation slowly and carefully. [Otherwise, there is a risk of device damage, unintended site damage, residual harvested tissue in the body, crush injury, or pain in patients.]
- Do not deflect the needles, sticks, or the Surelock during operation. [Otherwise, there is a risk of device damage.]
- Be careful not to let the needlepoint touch anywhere other than the intended contact surface. [Otherwise, there is a risk of the needlepoint damage.]
- Do not use the device to inject or drain fluids into or out of the body.
- After the procedure, take appropriate measures such as hemostasis.
- When wiping off the inner or outer needle, be careful not to get needle stick injuries and do not overload the device. Also, ensure that no foreign substance adheres to the needle parts. [Otherwise, there is a risk of device damage.]
- When performing repeated biopsies, check the appearance and ensure no damage (scratches, cracks, deformation, discoloration, and others) before use.


10.3 Other precautions

- The user must report all adverse events resulting from the use of the device to the manufacturer and to the competent authorities of the member state where the user or patient is.
- After use, dispose of the device safely and with care to prevent the risk of infection, following institutional, local, regional and national guidelines.


11 Device specifications

11.1 Configuration (typical diagram)


BMB body


BMB stylet


stick

surelock

safety cap

△ mark: Shown the cutting side of pipe.

11.2 Material

Raw material for parts that contact the blood, body fluids, mucous membranes, and others is below.

- Inner and outer needles, stick, surelock: stainless steel

The device is gamma-sterilized.

12 Instructions for use

12.1 Procedure of bone marrow puncture

- 1) Ensure that the BMB stylet is set on the BMB body.
- 2) Hold the outer needle handle and the inner needle knob tightly and puncture the needle into the bone marrow cavity. Be careful of over-puncturing, as resistance decreases when entering the bone marrow cavity.
- 3) After insertion into the bone marrow cavity, pull the BMB stylet straight out against the outer needle.


12.2 Procedure of bone marrow aspiration

If aspiration is necessary, perform it before the biopsy.

- 4) Connect a sterile syringe to the outer needle handle and perform bone marrow aspiration. Use a syringe that comply with ISO 80369-7.
- 5) Hold the outer needle handle when removing the syringe after aspiration is complete.
- 6) Remove the BMB body and hemostasis.

12.3 Procedure of bone marrow biopsy


- 7) Perform steps 1) to 3) for puncture.
- 8) Advance the BMB body while rotating it alternately from side to side to capture tissue into the outer needle.
- 9) Insert the stick into the BMB body and advance it to the point where it reaches the captured tissue. You can check the approximate length of the captured tissue at this point by measuring the scale of the marker.


- 10) After checking the length of the tissue, remove the stick.

12.4 Bone marrow sampling**With the Surelock**

- 11) Holding the outer needle handle, insert the Surelock completely. Do not rotate the Surelock during insertion. [Otherwise, there is a risk of device breakage.]
- 12) Rotate the BMB body 360 degrees to cut off the bone marrow sample. Rotate the BMB body, not the Surelock alone. [Otherwise, there is a risk of device breakage.]
- 13) Remove the outer needle from the puncture site while rotating the BMB body half-turn alternately from side to side. Do not pull the Surelock alone.
- 14) Confirm that the Δ mark on the hub part of the surelock is on the top, and pull out the surelock from BMB body.
- 15) Insert the stick into the Surelock from the side with the hub and push out the sample.

**Without the Surelock**

- 11) Turn the BMB body about ten times 360 degrees while performing a left to right motion.
- 12) After retracting the outer needle by approximately 2 to 3 mm, slightly change the puncture angle and advance the outer needle again by roughly 2 to 3 mm.
- 13) Turn the BMB body 360 degrees to cut off the sampling tissue and capture the sampling tissue in the outer needle.
- 14) Remove the outer needle from the puncture site while rotating the BMB body half-turn alternately from side to side.
- 15) Connect the safety cap with the larger opening to the tip of the outer needle.
- 16) Firmly insert the stick through the smaller opening of the safety cap and push the tissue out of the BMB body.

**12.5 At the end of the biopsy**

Reinsert the BMB stylet into the BMB body and connect the safety cap to the tip of the needle for safety.

13 Storage

- Keep dry. Keep away from direct sunlight, high temperature, and humidity. Store in a clean place.

14 Expiry date

- Refer to the expiry date printed on each package.
- Store properly and do not use the device after the expiry date.

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