Instructions for use**TSK STARCUT BIOPSY NEEDLE**
TSK STARCUT BIOPSY NEEDLE with Co-axial Guide Needle

The instructions for use covers TSK STARCUT BIOPSY NEEDLE (hereinafter referred to as "the device") and TSK STARCUT BIOPSY NEEDLE with Co-axial Guide Needle [with TSK Co-axial Guide Needle (hereinafter referred to as "introducer needle") attached].

1 Intended purpose

The device is a biopsy needle used to puncture human body and collect tissue for inspection, treatment, or diagnosis.

2 Indications

Patients requiring biopsy sample.

3 Intended users

The device is for use by medical doctors or health care professionals skilled in the procedures, which means it is not for use by anyone else. Only use by those medical doctors or health care professionals.

4 Target patients

The device is intended for use in patients 11 years of age and older of any gender, excluding cases of contraindications. The device should not be used in pregnant and lactating women.

5 Performance-characteristics

Insert into the human body to collect a range of soft tissue samples, excluding central nervous systems, for inspection, diagnosis, or treatment of patients.

6 Clinical benefits

The biopsy needle is a simple ancillary device that is widely used in many clinical procedures; the clinical benefit of the device is to have a positive impact on patient management:

- Allow management of patients' medical conditions by providing a means of collecting soft tissue samples locally in order to meet patients' medical needs.
- Allow medical professionals to manage medical conditions that require collection of soft tissue samples for diagnostic purposes.

7 Undesirable side effects

Possible adverse events associated with the use of the device are listed in the following:

- Air embolism
- Arteriovenous fistula due to puncture
- Bleeding
- Brain infarction and myocardial ischemia caused by aeroembolism
- Bruising
- Edema
- Hematoma
- Hemorrhage
- Hemothorax
- Iatrogenic displacement of the tumor cells
- Infection
- Needle stick injury
- Pain/discomfort
- Pericardial effusion
- Pleuritis
- Pseudoaneurysm formation
- Subcutaneous emphysema
- Swelling/mass sensation
- Vasovagal reflex

8 Contra-indications

- Do not use for patients who are receiving anticoagulant therapy, or have bleeding disorders or bleeding tendency.
- Do not use for the heart, the central circulatory system and the central nervous system.

9 Prohibitions

- Do not reuse. [Otherwise, there is a risk of infection.]
- Do not resterilize. [Re-sterilization does not assure device sterility and could damage or degrade the device.]

10 Warning**10.1 General precautions**

1.	Do not use the device under MRI. [Otherwise, there is a risk of MRI accidents.]
2.	Before using the device, read and follow the instructions for use.
3.	Do not use the device unless you are a medical doctor or a health care professional skilled in the procedure.
4.	If the patient or the device shows any abnormality, discontinue use immediately and take necessary measures.
5.	Check the device's appearance before use to ensure no damage (scratches, cracks, deformation, discoloration, and others).
6.	Do not use the device if the sterile bag damage (scratched, cracks, deformation, discoloration, and others) or the unintended opening before use are shown.
7.	Refer to the instructions if you use the device with pharmaceuticals or other medical devices.
8.	The device could be contraindicated in patients receiving anticoagulant


	treatment or in patients with bleeding disorders or bleeding tendency.
9.	Patient management after biopsy depends on the biopsy technique and the physiological state of the individual patient. Monitor the patient's vital signs and strive to prevent and treat adverse events associated with the biopsy procedure.
10.	Do not use the device for any purpose other than those described in the "instructions for use."
11.	Do not modify the device.
12.	Operation practicing results in damage to the needlepoint and poor sharpness. Use samples for operation practice.

10.2 Precautions for use


- Open and remove the device with the appropriate aseptic technique. [Otherwise, there is a risk of infection.]
- Use the device immediately after opening the package. [Otherwise, bacteria could multiply and cause infection.]
- Check the selection of gauge size, length, and the stroke length of 11mm or 22mm in advance. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- Be careful when removing or attaching the protector. [Otherwise, there is a risk of needle stick injury.]
- Be careful that the cutting edge does not contact the protector, biopsy needle guide, or other unintended contact surfaces. [Otherwise, there is a risk of the needlepoint damage.]
- Do not touch the needle directly with your hands. [Otherwise, there is a risk of infection.]
- The location of the puncture site should be confirmed by diagnostic equipment and palpation.
- Check the puncture route carefully when inserting the needle. [Otherwise, there is a risk of damaging unintended areas.]
- Do not deflect the needle during operation. [Otherwise, there is a risk of device damage.]
- When using an introducer needle, use an introducer needle that is compatible with the device. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- When using a biopsy needle guide such as an ultrasonic probe in combination, do not forcibly set the needle on a guide whose shape does not fit it. [Otherwise, there is a risk of device damage.]
- Do not use the device to inject or drain fluids into or out of the body.
- If the needle does not advance as desired during insertion, do not forcibly advance it. [Otherwise, there is a risk of device damage.]
- The amount of tissue sampling depends on the nature of the target tissue.
- There have been reports of air embolism due to puncture operations in the lungs, leading to cerebral infarction and myocardial ischemia. Thus, if the patient develops any of these symptoms, keep the patient's head in a low position, diagnose with CT and others, and take appropriate measures immediately. In severe cases, consider hyperbaric oxygen therapy without delay.
- After the procedure, take appropriate measures such as hemostasis.
- When wiping off the inner or outer needle, be careful not to get needle stick injuries and do not overload the device. Also, ensure that no foreign substance adheres to the needle parts. [Otherwise, there is a risk of device damage.]
- When performing repeated biopsies, check the appearance and ensure no damage (scratches, cracks, deformation, discoloration, and others) before use.

10.3 Other precautions

- The user must report all adverse events resulting from the use of the device to the manufacturer and to the competent authorities of the member state where the user or patient is.
- After use, dispose of the device safely and with care to prevent the risk of infection, following institutional, local, regional and national guidelines.

11 Device specifications**11.1 Configuration (typical diagram)****TSK STARCUT BIOPSY NEEDLE****TSK Co-axial Guide Needle**

Only model of TSK STARCUT BIOPSY NEEDLE with Co-axial Guide Needle is attached.

11.2 Material

Raw material for parts that contact the blood, body fluids, mucous membranes, and others is below.

- Inner and outer needles: stainless steel

The device is gamma-sterilized.


12 Instructions for use**12.1 Procedure of using introducer needle**

If not use the introducer needle, perform from 4).

- 1) Remove the protector of introducer needle.
- 2) Puncture the introducer needle towards the marginal zone of the target site. If there is resistance and the needle does not penetrate as expected during insertion, do not forcibly insert the needle. [Otherwise, there is a risk of device damage.]
- 3) Remove the inner needle from the outer needle of the introducer needle. When removing the inner needle, hold the outer needle hub and turn the inner needle hub to unlock it.

12.2 Stroke length selection

- 4) Remove the safety. Then hold the TSK hub and remove the protector. Do not remove the protector first. [Otherwise, there is a risk of needle stick injury.]
 - 5) Pull the Trigger until locked to store stress in the spring. There are two levels of locking. The locking position of the Trigger selects the stroke length of the inner needle during sample collection. A window on the body shows the locking position of the Trigger. When pulling the Trigger, do not apply any load such as bending. [Otherwise, there is a risk that the inner needle could come off and damage the device.]
- At the first lock, the Trigger locks in position "10," and the stroke length of the inner needle during sample collection is 11 mm.
 - At the second lock, the Trigger locks in position "20," and the stroke length of the inner needle during sample collection is 22 mm.


Puncture with the Trigger pulled.


12.3 Puncture

If not use the introducer needle, perform from 7).


- 6) Insert the needle of the device into the outer needle of the introducer needle.
- 7) With the Trigger pulled, puncture to the marginal zone of the biopsy site. Do not push the Trigger until the needle reaches the intended position. [Otherwise, there is a risk of damaging the unintended site.]

**12.4 Tissue sampling**

- 8) Push the Trigger to expose the inner needle.
- 9) When you push the lower circular part of the Trigger (where "Push here" is printed) to the end, the outer needle automatically fires, and the tissue-sampling groove samples the tissue.

**12.5 Tissue removal**

- 10) (i) After needle removal, pull Trigger back until Trigger is locked in the "10" position; do not perform tissue removal operation with Trigger in the "10" position. [Otherwise, there is a risk of needle stick injury.]
- 11) (ii) Then push the Trigger to expose the tissue-sampling groove and remove the tissue. Do not over-push the Trigger and cause the outer needle to fire. [Otherwise, there is a risk of needle stick injury.]

- 12) If you repeat the biopsy, return to 12.2 Stroke length selection.

12.6 Removal the introducer needle

If not use the introducer needle, the following procedure is not performed.

- 13) Insert the inner needle into the outer needle of the introducer needle. Do not remove the outer needle alone. [Otherwise, there is a risk of device damage.]
- 14) While holding the outer needle hub, turn and lock the inner needle hub.
- 15) Remove the introducer needle.

13 Storage

- Keep dry. Keep away from direct sunlight, high temperature, and humidity. Store in a clean place.

14 Expiry date

- For the useful life of the device that has been stored properly as per the storage instructions, refer to the use-by date printed on each individual package.
- Store properly and do not use the device after the expiry date.

15 Manufacturer's information

TSK Laboratory, Japan
1510-1 Soja-Machi, Tochigi-Shi, Tochigi-Ken, 328-0002 JAPAN
TEL : +81-282-27-0005
FAX : +81-282-25-6511
URL : www.tsklab.co.jp/en/
E-mail : tsk@tsklab.co.jp

Anwendungshinweise**TSK STARCUT BIOPSY NEEDLE**
TSK STARCUT BIOPSY NEEDLE with Co-axial Guide Needle

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