Instructions for use

TSK Co-axial Guide Needle**1 Intended purpose**

The TSK Co-axial Guide Needle is a single use puncture device, intended to assist the placement of a procedural needle into the body for treatment and/or diagnosis purposes.

2 Indication

The device is indicated for patients requiring biopsy sample.

3 Intended users

The device is intended for use by doctors or medical professionals skilled in the procedure.

4 Target patients

The device is intended for use in patients 11 years of age and older of any gender, excluding cases of contraindications. The device should not be used in pregnant and lactating women.

5 Performance characteristics

When used together with TSK biopsy needles the device enables introduction of the procedure needle into the patient's body. The device provides an unobstructed channel for collecting tissue samples from a patient through the outer needle lumen of the TSK biopsy needle.

6 Clinical benefits

By virtue of enabling use of TSK biopsy needles, the device provides an indirect clinical benefit.

The indirect clinical benefit of the device is to have a positive impact on the management of the patient's medical condition by providing an unobstructed channel for collecting tissue samples.

7 Undesirable side effects

Possible adverse events associated with the use of the device are listed in the following:

- pneumothorax
- hemoptysis
- hemothorax
- air embolism
- ascites / hemoperitoneum
- bleeding / hemorrhage / hematoma
- perforation
- seeding of cancer cells along needle tract
- pain

8 Contraindications

- Do not use for patients who are receiving anticoagulant therapy, or have bleeding disorders or bleeding tendency.
- Do not use for the heart, the central circulatory system and the central nervous system.

9 Prohibitions

- Do not reuse. [Otherwise, there is a risk of infection.]
- Do not resterilize. [Re-sterilization does not assure device sterility and could damage or degrade the device.]

10 Warning**10.1 General precautions**

1.	Do not use the device under MRI. [Otherwise, there is a risk of MRI accidents.]
2.	Before using the device, read and follow the instructions for use.
3.	Do not use the device unless you are a medical doctor or a health care professional skilled in the procedure.
4.	If the patient or the device shows any abnormality, discontinue use immediately and take necessary measures.
5.	Check the device's appearance before use to ensure no damage (scratches, cracks, deformation, discoloration, and others).
6.	Do not use the device if the sterile bag damage (scratched, cracks, deformation, discoloration, and others) or the unintended opening before use are shown.
7.	Refer to the instructions if you use the device with pharmaceuticals or other medical devices.
8.	Patient management after the procedure depends on technical details and the physiological state of the individual patient. Monitor the patient's vital signs and strive to prevent and treat adverse events associated with the procedure.
9.	Do not use the device for any purpose other than those described in the "instructions for use."
10.	Do not modify the device.


- | | |
|-----|---|
| 11. | Operation practicing results in damage to the needlepoint and poor sharpness. Use samples for operation practice. |
|-----|---|

10.2 Precautions for use

- Open and remove the device with the appropriate aseptic handling. [Otherwise, the device could become contaminated and cause infection.]
- Use the device immediately after opening the package. [Otherwise, bacteria could multiply and cause infection.]
- Protect the hub of the device from ethanol and organic solvents. [Otherwise, there is a risk of hub damage.]
- Confirm needle gauge size and needle length selection in advance. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- Before use, ensure that the medical device used in conjunction with the device can pass through the outer needle. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- The location of the puncture site should be confirmed by diagnostic equipment and palpation.
- Check the puncture route carefully when inserting. [Otherwise, there is a risk of injuring an unintended site.]
- Do not touch the inner or outer needle directly with your hand. [Otherwise, there is a risk of infection.]
- The device is not intended for use in stand-alone operations or procedures.
- Do not use the device to inject fluids into or drain fluids from the body.
- Be careful when removing or attaching the protector. [Otherwise, there is a risk of needle stick injury.]
- Be careful that the needlepoint does not contact unintended contact surfaces. [Otherwise, there is a risk of needlepoint damage.]
- Ensure that the inner and outer needle hubs are correctly set. Hold these hubs all together during puncture. [Otherwise, the inner needle could be dislodged, failing to puncture as expected.]
- Do not deflect the needle during puncture, removal, and others. [Otherwise, there is a risk of device damage.]
- Do not pull back the device used in combination while in the outer needle after insertion of the combination device (guidewire, catheter, and others). [Otherwise, there is a risk of combination device damage.]
- After the procedure, take appropriate measures such as hemostasis.

10.3 Other precautions

- The user must report all adverse events resulting from the use of the device to the manufacturer and to the competent authorities of the member state where the user or patient is.
- After use, dispose of the device safely and with care to prevent the risk of infection, following institutional, local, regional and national guidelines.

11 Device specifications**11.1 Configuration (typical diagram)****11.2 Material**

Raw material for parts that contact the blood, body fluids, mucous membranes, and others is below.

- Inner and outer needles: stainless steel

The device is ethylene oxide sterilized.

12 Instructions for use

- 1) Disinfect the puncture site.
- 2) Ensure that there are no scratches or bends on the device before use.
- 3) Remove the protector from the needle of the product.
- 4) Provide appropriate local anesthesia.
- 5) Puncture towards the marginal zone of the biopsy site. If there is resistance and the needle does not advance as expected during insertion, do not forcibly insert the needle. [Otherwise, there is a risk of device damage.]
- 6) Remove the inner needle of the device from the outer one, and perform the intended procedure. When removing, turn the "inner needle hub" to unlock.
- 7) After completion of the procedure, remove the outer needle carefully.

13 Storage

- Keep dry. Keep away from direct sunlight, high temperature, and humidity. Store in a clean place.

14 Expiry date

- For the useful life of the device that has been stored properly as per the storage instructions, refer to the use-by date printed on each individual package.
- Store properly and do not use the device after the expiry date.

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