Instructions for use
TSK Syringe

1 Intended purpose

The device is a general purpose syringe intended for use with a needle to inject or withdraw substances. The device is also indicated for intravitreal use.

2 Indication

The device is intended for use with patients that need injection or withdrawal of substances for examination, treatment, diagnosis, or prevention.

3 Intended users

The device is intended for use by healthcare professionals skilled in the procedure.

4 Target patients

The device is intended to be used in pediatric and adult populations (including pregnant and breastfeeding women), regardless of gender, and who do not meet contraindications.

5 Performance characteristics

- Connection with a needle
- Provide measures for withdrawal and injection of substances

6 Clinical benefits

- Allow management of patients' medical conditions by providing a means of delivering essential substances locally or systemically in order to meet patients' medical needs.
- Allow medical professionals to manage medical conditions that require removal of biological fluid, and obtain biological fluid samples for diagnostic purposes.

7 Undesirable side effects

Possible adverse events associated with the use of the device are listed in the following:

- hematoma
- infection
- inflammation, swelling
- bleeding, internal bleeding
- vascular embolism
- ocular inflammation
- floaters

8 Contraindications

The device can be used in any standard application for a syringe. There are no contraindications.

9 Prohibitions

- Do not reuse. [Otherwise, there is a risk of infection.]
- Do not resterilize. [Re-sterilization does not assure device sterility and could damage or degrade the device.]

10 Warnings

10.1 General precautions	
1.	Before using the device, read and follow the instructions for use.
2.	Do not use the device unless you are a healthcare professional skilled in procedure.
3.	If the patient or the device shows any abnormality, discontinue use immediately and take necessary measures.
4.	Check the device's appearance before use to ensure no damage (scratches, cracks, deformation, discoloration, and others).
5.	Do not use the device if the sterile bag is opened or damaged (scratched, cracks, deformation, discoloration, and others).
6.	Refer to the associated instructions if you use the device with pharmaceuticals or other medical devices.
7.	Do not use the device for any purpose other than those described in the "instructions for use."

8.	Do not modify the device.
9.	Use only with devices that are ISO 80369-7 compliant.

10.2 Precautions for use


- Use the device immediately after opening the package. [Otherwise, bacteria could multiply and cause infection.]
- Open the bag from the designated opening position.
- Do not touch the tip and inside of barrel with your hands. [Otherwise, there is a risk of infection.]
- Check the nominal capacity before use. [Otherwise, there is a risk that you could not perform the appropriate procedure.]
- Before use with a needle, ensure that the connection to the needle is not loose. [Otherwise, there is a risk of leakage or device dropout.]
- Do not apply excessive force when connecting to a needle. [Otherwise, there is a risk of device damage.]
- Ensure that there are no damage, loose connections, or leaks during use.
- Ensure that pharmaceuticals, blood, and other liquids do not enter the connection. [Otherwise, the connecting part could loosen.]
- When penetrating vial closures, do not apply lateral force on the device. [Otherwise, there is a risk of device damage.]
- Do not push the plunger with excessive force during the procedure. [Otherwise, there is a risk of injury to the patient.]
- When aspirating, do not pull the plunger diagonally. [Otherwise, there are risk that air gets mixed in.]
- When aspirating, do not pull the plunger too far so that it comes out the back of the barrel.
- When injecting substances into patients, fill the entire needle tube with the substances before use.
- Protect the hub of the device from ethanol and organic solvents. [Otherwise, there is a risk of device damage.]
- Do not fill with substances that have a temperature significantly different from the ambient environment. [Otherwise, there are risks that the filling volume is incorrect, and/or operability is deteriorated due to material shrinkage.]
- Do not rub the printed markings excessively. [Otherwise, there is a risk that readability will deteriorate.]
- Do not use for high-pressure injection, such as of a contrast medium. [Otherwise, there are risks including leakage and damage.]

10.3 Other precautions

- The user must report all adverse events resulting from the use of the device to the manufacturer and to the competent authorities of the member state where the user or patient is.
- After use, dispose of the device safely and with care to prevent the risk of infection, following institutional, local, regional and national guidelines.
- Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

11 Device specifications

11.1 Configuration (typical diagram)

The device complies with ISO 7886-1 and ISO 80369-7.

11.2 Materials

Raw materials for parts that contact the blood, body fluids, mucous membranes, and other patient tissues are:

- [1] barrel: copolyester resin
- [2] gasket: compound resin
- lubricant: silicone

The device is gamma-sterilized.

12 Instructions for use

- 1) If necessary, wear gloves in advance.
- 2) Open the sterile bag using appropriate aseptic technique and remove the device.
- 3) Connect firmly to the needle or other device.
- 4) When filling with substances, the substances can be aspirated by pulling the plunger.
- 5) Disinfect the puncture site and perform the puncture.
- 6) Perform the injection or withdrawal of substances.
- 7) Remove the needle and/or other devices, and if necessary, stop the bleeding.
- 8) Dispose of the device following standard safety protocols.

13 Storage

- Keep dry. Keep away from direct sunlight, high temperature, and humidity. Store in a clean place.

14 Expiry date

- For the useful life of the device that has been stored properly as per the storage instructions, refer to the use-by date printed on each individual package.
- Store properly and do not use the device after the expiry date.

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