



TSK

Product Catalog

TSK OptiXion™

Better precision. Sharper vision.


(1) TSK OptiXion™ NEEDLE

HUB MARKING

to identify the direction of the bevel in situations where bevel orientation is particularly hard to see.


WIDE LUMEN

allows smoother injection.





TSK OptiXion™ Needle is a single-use injection needle mainly used to administrate medications and to drain intraocular substances. The tip of **TSK OptiXion™ Needle** is designed to reduce penetration force when punctured into the eye.

(2) TSK OptiXion™ Viscoelastic Cannula



TSK OptiXion™ Viscoelastic Cannula is a single-use injection cannula mainly used for inject viscoelastic substances (OVD) into the eye.

Product Lineup

[PRODUCT TYPE]	[TIP SHAPE]	[SIZE]	[PRODUCT CODE]
(1) TSK OptiXion™ NEEDLE		● 34G x 9mm^{ETW} ● 33G x 9mm^{ETW} ● 32G x 9mm^{ETW} ● 32G x 13mm^{ETW} ● 30G x 13mm^{ETW}	OTX-018009E OTX-020009E OTX-023009E OTX-023013E OTX-030013E
(2) TSK OptiXion™ Viscoelastic Cannula		● 27G x 22mm ● 25G x 22mm	VSC-040022 VSC-050022

※ Please refer to the latest package insert before use.

※ Some products may not be available for sale in certain countries.
For further information, please inquire your local distributor.

50pcs/ box
Gamma ray sterilization, single-use

Distributed by:

TSK TSK LABORATORY, JAPAN

1510-1 Soja-machi, Tochigi-shi, Tochigi-ken, 328-0002 Japan
TEL: +81 282 27 0005 FAX: +81 282 25 6511

www.tsklab.co.jp/en/