

TSK

TSK SELECT™ NEEDLE

TSK
HYPODERMIC
NEEDLE





Needle hub marking for easier identification

The **TSK SELECT™ Needle** features a hub with marking that helps practitioner to identify the direction of the cutting edge. In the case where it is difficult to see the cutting edge due to the ultra-fine 34G needle, or when the tip of the needle has penetrated the skin, the direction of the cutting edge can be easily confirmed for precise and safe handling.



34G (OD: 0.18mm)

34G ultra-fine needle reduces burden on patients

TSK SELECT™ Needle's thinnest gauge, the "true*" 34G (outer diameter of 0.18mm), can be used in areas where discomfort during puncture should be reduced, and at times when recovery after a procedure is of particular concern. Furthermore, the tip of **TSK SELECT™ Needle** is designed to maintain sharpness even after multiple punctures.

*Compliant with International Organization for Standardization (ISO) standards.

Universal hub design: works with both slip type and luer-lock type syringes*

When paired together with lock-type syringes, such as the TSK Syringe, it can further reduce the risk of pop-off and leakage that can occur when injecting high-viscosity solutions or when using thin needle.

*Only syringes compliant with International Organization for Standardization (ISO) standards.



Low dead space design minimizes residual chemical solution

The hub of **TSK SELECT™ Needle** is designed to minimize residual solution left in the needle hub after injection, maximizing the use of solutions. When paired with **TSK Syringes**, it can further lessen the dead space inside the hub.



More information
on TSK Syringe



TSK SELECT™ NEEDLE | Product Lineup

PRODUCT CODE	PRODUCT NAME
TSN-018004	TSK SELECT NEEDLE 34G (0.18mm) x 4mm
TSN-018009	TSK SELECT NEEDLE 34G (0.18mm) x 9mm
TSN-020004R	TSK SELECT NEEDLE 33G (0.20mm) x 4mm RW
TSN-020009R	TSK SELECT NEEDLE 33G (0.20mm) x 9mm RW
TSN-020004	TSK SELECT NEEDLE 33G (0.20mm) x 4mm
TSN-020009	TSK SELECT NEEDLE 33G (0.20mm) x 9mm
TSN-020013	TSK SELECT NEEDLE 33G (0.20mm) x 13mm
TSN-023004	TSK SELECT NEEDLE 32G (0.23mm) x 4mm
TSN-023009	TSK SELECT NEEDLE 32G (0.23mm) x 9mm
TSN-023013	TSK SELECT NEEDLE 32G (0.23mm) x 13mm
TSN-025004	TSK SELECT NEEDLE 31G (0.25mm) x 4mm
TSN-025009	TSK SELECT NEEDLE 31G (0.25mm) x 9mm
TSN-025013	TSK SELECT NEEDLE 31G (0.25mm) x 13mm
TSN-030004	TSK SELECT NEEDLE 30G (0.30mm) x 4mm
TSN-030009	TSK SELECT NEEDLE 30G (0.30mm) x 9mm
TSN-030013	TSK SELECT NEEDLE 30G (0.30mm) x 13mm
TSN-030013E	TSK SELECT NEEDLE 30G (0.30mm) x 13mm ETW
TSN-030025	TSK SELECT NEEDLE 30G (0.30mm) x 25mm
TSN-040013	TSK SELECT NEEDLE 27G (0.40mm) x 13mm
TSN-040019R	TSK SELECT NEEDLE 27G (0.40mm) x 19mm RW
TSN-050013	TSK SELECT NEEDLE 25G (0.50mm) x 13mm
TSN-050025R	TSK SELECT NEEDLE 25G (0.50mm) x 25mm RW

※Please refer to the latest package insert before use.

Quantity: 100pcs/box

Trade Name: TSK SELECT™ NEEDLE

Product Name: TSK HYPODERMIC NEEDLE

Gamma ray sterilization , Single-use

※Some products may not be available for sale in certain countries.
For further information, please inquire your local distributor.

Distributed by:



TSK LABORATORY, JAPAN

1510-1, Soja-Machi, Tochigi-Shi, Tochigi-Ken
328-0002 JAPAN
ITL TEL +81 282 27 0005
ITL FAX +81 282 25 6511
CAT-TSN-E03 2025.7